

Mimea 6th Generation

THE NEW MEASUREMENT AND CONTROL SYSTEM

KEY FEATURES

- Advanced high-speed data processing architecture
- A new standard in signal quality for data acquisition and signal generation
- Unparalleled flexibility in signal modulation and mixing
- Optimal I/O placement capability via galvanically isolated optical interconnects
- Ready and open for AI integration



Nanonis™

A member of SPECSGROUP

Mimea 6th Generation

A QUANTUM LEAP FOR ADVANCED EXPERIMENTS

Made for scientific excellence

The sixth-generation Nanonis Mimea SPM control system represents a new benchmark for scanning probe microscopy, built upon a legacy of over 20 years of intensive research and development and the experience gained from over 1600 systems delivered to the global research community.

For more than a decade, the fifth generation Nanonis Mimea SPM control system has been the primary choice for the world's most demanding SPM experiments, serving as the technical foundation for countless scientific breakthroughs and high-impact publications. Following our core philosophy that the measurement instruments should never be the limiting factor in a scientific experiment, the system has undergone continuous product development driven directly by the needs of the SPM community. This unique collaboration has resulted in over a hundred new features and specialized solutions, developed by an in-house software engineering team determined to provide effective solutions and new features in a short time.

However, as experimental requirements push the boundaries of physics, the hardware platform of the controller must follow suit. The sixth generation Nanonis Mimea SPM control system addresses current and future demands through a ground-up redesign. And while the external appearance might seem familiar, the completely new hardware design and software architecture are capable of delivering performance levels previously considered inconceivable. Paired with a higher degree of modularity, the new architecture lays a solid foundation for further decades of continuous development. This major investment ensures that the system is not only capable of standard SPM measurements but is also flexible enough to address emerging applications and future scientific challenges, providing greater value for money.

By combining advanced technology and community-driven development, the Nanonis Mimea SPM control system continues to be a major contributor to global scientific excellence.

Building on legacy, developing for the future

- 1600 systems installed worldwide
- 25 years of continuous research & development

Modularity and serviceability

- Compatible with existing high voltage electronics and adaptation kits
- Modular system and modular instruments: easy to upgrade and to service
- Extensive internal monitoring to ensure smooth operation

Ultrastable, ultra low-noise and fully monitored signal generation

- Voltage or current source mode
- 31-bit DAC architecture
- 6× higher maximum bandwidth
- Up to 6× lower output noise

APIs, scripting and STM simulator

- Ready for AI: platform-agnostic APIs and native Python interface

A state of the art 25 MHz AFM interface

- Arithmetic unit for full frequency flexibility
- 2× the number of modulators and demodulators
- 5× larger bandwidth
- 5× lower input noise
- Up to 30× lower output noise

State of the art data acquisition in any experimental situation

- Selectable input gain of 1, 10 or 100
- Selectable AC coupling
- Up to 20× less input noise
- Up to 8× higher input resolution
- 100× higher input impedance
- 2× higher sampling rate

Designed for high speed and big data

- 2× larger number of physical and software signals
- Up to 16 lock-in demodulators
- 6× higher data transfer rate for each instrument
- 5× higher sampling rates possible
- Up to 50m galvanically isolated optical interconnects between instruments





RC6 Real-Time Controller

A key benefit of the Nanonis high-speed modular architecture is the direct access and visualization of all raw signals available in the system, eliminating the need for pre-processing and data-rate reduction in add-on instruments. To ensure a low-latency and reliable data transfer at a significantly higher speed than before, a revolutionary data transfer architecture had to be designed. It supports a 7× higher data transfer rate, sampling rates up to 200 MS/s, synchronous data input and output over all interfaces in the system, and plug&play operation.

Designed for distributed measurements and optimized for low HF signature and small signals, the architecture enables a placement of the signal interfaces at a significant distance from the central processing core, the RC6, using galvanically isolated optical fibers. For complex measurements, a single system can handle over 128 outputs or inputs and over 250 monitoring signals, with over 200 signals simultaneously available in the user interface. To handle this enormous amount of data, the RC6 offers increased processing power with a larger FPGA, a more powerful CPU and faster storage.

SC6 Signal Conversion

The SC6 redefines the meaning of signal interface, offering state-of-the-art performance in any experimental situation. The 24-bit 2 MS/s ADCs paired with a low-noise, low-drift input stage with variable gain and 200 kHz bandwidth, offer a 3× reduction in broadband noise and an 8× higher effective resolution compared to the previous generation, with oversampling effective down to 1 Hz. A high-impedance input buffer, true-differential inputs and flexible GND schemes ensure optimal operation with any sample and electrical configuration, while selectable AC-coupling with high gain allows the detection of smallest modulations.

The 31-bit DAC architecture delivers incredibly clean signals down to smallest amplitudes eliminating the need for external adders. Variable filters offer a 6× higher bandwidth and an optimized impulse response for stable feedback. Output noise is twice as low, or even 6× lower in the $\pm 2.5\text{V}$ output range, with the improvement carrying over to frequencies below 0.1 Hz. Output current and voltage are constantly monitored, and all outputs can act as regulated current sources with maximum currents up to 25 mA.



OC6 Oscillation Controller

The OC6 is a state-of-the-art signal generator with 5× the bandwidth compared to the OC4, capable of generating small amplitude signals up to 100MHz. The redesigned output signal path offers improved linearity, signal purity, and fully digital amplitude control. With user-selectable ranges and filters to reduce broadband noise and optimize signal quality for the application, it offers uncompromised performance also for signals down to DC. The high-current output stage can be bypassed for applications not requiring large driving currents, lowering noise to less than 5nV/sqrt(Hz).

The OC6 is also a sophisticated digitizer, with a DC-coupled high-impedance, fully differential input stage, able to handle signals from DC to 25MHz or even 100MHz with highest accuracy. Great care has been taken to avoid any interference with even the smallest input signals, resulting in 5× less noise than the predecessor, a significantly lower 1/f corner frequency and lower distortion. The OC6 provides researchers with a unique, all-in-one tool designed to meet the most demanding requirements in AFM measurements, but able to go significantly beyond those.

Software and add-on modules

The Nanonis V6 software is both familiar and significantly more powerful. It builds on the proven and reliable V5e software, but supports all the new hardware functionality, including the extensive monitoring and safety features and the over 200 signals accessible in real-time. It offers an unparalleled flexibility in signal modulation and mixing, and twice the number of frequency generators and demodulators.

It integrates specialized tools derived from the Nanonis Tremea quantum transport measurement platform. Together with embedded control of external instruments like magnet and temperature controllers or RF generators, such modules make combined SPM and transport measurements much more manageable.

The Programming Interface, a game-changing tool for researchers for over two decades, now runs on a generic framework with dedicated Python and LabVIEW interfaces for immediate productivity, and is ready for any advances in machine-learning based routines. The V6 software combines reliability with modular features and continuous development.

Technical Data

RC6 Real-Time Controller

Computing	6-core CPU, 8GB RAM, 512GB SSD, FPGA-card
Operating system	Linux RT
FPGA card	PXIE, high-speed serial interface
Connectivity	4× SC6 max., 2× OC6 dual max.

SC6 Signal Conversion

Power supply	Powered from RC6
Electrical GND	Floating up to 30V or PE, remote controlled, ground loop detection

Input

Hardware interface	8× BNC connectors, differential over 2 connectors, pseudodifferential, single-ended Shield on PE, GND or floating, user-configurable
Differential input range	±10V, ±1V, ±100mV (gain 1, 10, 100)
Input impedance	100MΩ, > 1GΩ, user configurable
Input capacitance	<40pF
Input bias current	1.5pA typ., 2pA max.
Coupling	DC or AC, user-configurable
Analog bandwidth	DC-200kHz (-3dB) 6 th order Bessel low-pass filter
AD-converter	24-bit 2MS/s
Effective resolution	>25 bit with >4000× oversampling
INL / DNL	±2LSB / ±0.5LSB typical
Input noise density at gain 100	<7nV/√Hz > 100Hz <20nV/√Hz @ 1Hz
Measurement noise at gain 100	<5μVrms @ 2MS/s <0.5μVrms @ 15kS/s <80nVrms @ 120S/s
0.1–10Hz noise at gain 100	<60nVrms
SFDR, 50% amplitude signal, 1 kHz	>125dB @ gain 1 >105dB @ gain 10 >105dB @ gain 100

Output

Hardware interface	8× BNC connectors
Output voltage ranges	±10V, ±2.5V
Output current ranges	±25mA, ±12.5μA
Output impedance	<1Ω (25mA), 1kΩ (12.5μA)

Output

Analog bandwidth	DC-250kHz (-3dB), user-configurable 100kHz, 20kHz, 350Hz low-pass filters
DA converter	31-bit DAC architecture with active glitch compensation, 1MS/s
Effective resolution	23-bit typ., >24-bit with 350Hz filter at full sampling rate
INL / DNL	<±2ppm max. / <±1ppm max.
Output noise density ±10V range	<11nV/√Hz > 1kHz <25nV/√Hz @ 1Hz
Output noise density ±2.5V range	<3.5nV/√Hz > 1kHz <10nV/√Hz @ 1Hz
Output noise 10 – 300 kHz ±10V range (250 kHz / 100 kHz / 20 kHz / 350 Hz filter)	<6μVrms / <4.5μVrms / <2.5μVrms / <0.3μVrms
Output noise 0.1 Hz – 10 Hz ±10V range	<180nVrms
Output noise 0.1 Hz – 10 Hz ±2.5V range	<50nVrms <50nVrms
12h-Drift	<25μV p-p (@ 9.9V)

OC6 Oscillation Controller

Power supply	Powered from RC6
Electrical GND	Floating up to 30V or PE, remote controlled, ground loop detection

Input

Hardware interface	2× BNC connectors, differential or single-ended
Differential input range	±10V, ±3.3V, ..., ±1mV (gain 0.1, 0.3, ..., 1000)
Input impedance	50Ω, 100MΩ, > 1GΩ, DC-coupling 50Ω, 10MΩ, AC-coupling user configurable
Input capacitance (single-ended)	15pF (low-capacitance mode) 40pF (low-noise mode)
Input bias current (single-ended, input A, 60 °C)	<35pA typ. (low-capacitance mode) <140pA typ. (low-noise mode)
Coupling	DC or AC, user-configurable
Analog bandwidth	DC-25MHz, -5MHz, -500kHz (-3dB), 3 rd order Bessel low-pass filters, user selectable
AD-converter	16-bit 200MS/s
Effective resolution	>17bit with >256× oversampling, 1 – 10V range

Input	
INL / DNL	$\pm 1.6\text{LSB} / \pm 0.4\text{LSB}$ typical
Input noise density in 10mV range or lower, low-noise, single-ended	$< 3.5\text{nV}/\sqrt{\text{Hz}}$ > 100kHz $< 15\text{nV}/\sqrt{\text{Hz}}$ @ 100Hz
Measurement noise up to 500kHz, low-noise, single-ended, 128× oversampling	$< 22\mu\text{Vrms}$, 1V range $< 2.5\mu\text{Vrms}$, 10mV range or lower
Measurement noise up to 100MHz, low-noise, single-ended, no oversampling	$< 170\mu\text{Vrms}$, 1V range $< 28\mu\text{Vrms}$, 10mV range or lower
SFDR, 50% amplitude signal, 1 kHz	$> 70\text{dB}$, 10V range $> 90\text{dB}$, 1V range
Output	
Hardware interface	1 × BNC connector
Output voltage ranges	$\pm 10\text{V}$, $\pm 1\text{V}$, $\pm 100\text{mV}$, $\pm 10\text{mV}$
Output current	$\pm 200\text{mA}$ max. (high current) $\pm 100\text{mA}$ max. (low noise)
Output impedance	50Ω
Analog bandwidth	DC-25 MHz, -5 MHz, -500kHz (-3dB), 3 rd order Bessel low-pass filters, user selectable
DA converter	16-bit DAC, 200MS/s
INL / DNL	$< \pm 1\text{LSB} / \pm 0.5\text{LSB}$ typical
Output noise density $\pm 10\text{V}$ range, high current	$< 60\text{nV}/\sqrt{\text{Hz}}$ > 10kHz $< 800\text{nV}/\sqrt{\text{Hz}}$ @ 10Hz
Output noise density $\pm 1\text{V}$ range, high current	$< 25\text{nV}/\sqrt{\text{Hz}}$ > 10kHz $< 250\text{nV}/\sqrt{\text{Hz}}$ @ 10Hz
Output noise density $\pm 1\text{V}$ range, low noise	$< 7\text{nV}/\sqrt{\text{Hz}}$ > 10kHz $< 70\text{nV}/\sqrt{\text{Hz}}$ @ 10Hz
Output noise density $\pm 10\text{mV}$ range, low noise	$< 5\text{nV}/\sqrt{\text{Hz}}$ > 10kHz $< 25\text{nV}/\sqrt{\text{Hz}}$ @ 10Hz
Output noise up to 500kHz, low-noise (high-current)	$< 55\mu\text{Vrms}$ ($< 65\mu\text{Vrms}$), $\pm 10\text{V}$ range $< 6\mu\text{Vrms}$ ($< 15\mu\text{Vrms}$), $\pm 1\text{V}$ range $< 3.5\mu\text{Vrms}$ ($< 14\mu\text{Vrms}$), $\pm 10\text{mV}$ range
Output noise up to 100MHz, low-noise (high-current)	$< 55\mu\text{Vrms}$ ($< 240\mu\text{Vrms}$), $\pm 10\text{V}$ range $< 55\mu\text{Vrms}$ ($< 65\mu\text{Vrms}$), $\pm 1\text{V}$ range $< 55\mu\text{Vrms}$ ($< 62\mu\text{Vrms}$), $\pm 10\text{mV}$ range
Output noise 0.1 Hz – 10 Hz, low-noise (high current)	$< 5\mu\text{Vrms}$ ($< 6.5\mu\text{Vrms}$), $\pm 10\text{V}$ range $< 0.6\mu\text{Vrms}$ ($< 1.8\mu\text{Vrms}$), $\pm 1\text{V}$ range $< 0.2\mu\text{Vrms}$ ($< 1.5\mu\text{Vrms}$), $\pm 10\text{mV}$ range

Digital lines	
Ports	2 × 8 bidirectional lines on two D-sub 9 female connectors for each signal interface
Signal	3.3V TTL, max. 25mA per line
Maximum sampling frequency	250 kHz – 1 MHz depending on configuration
High speed digital lines	
Ports	2 × inputs and 2 × outputs on SMB male connectors for each signal interface
Signal	3.3V TTL, max. 24mA per line
Maximum sampling frequency	200MHz
Clock	
Ports	1 × input, 1 × output for active clock source
Frequency	10MHz, 100MHz input, 10MHz, 100MHz, 1GHz output square wave, 3.3V
Frequency calibration	$\pm 25\text{ppm}$ (standard VCXO) $\pm 1\text{ppm}$ (optional OCXO)
Frequency stability vs. temperature	$\pm 15\text{ppm}$ (standard VCXO) $\pm 20\text{ppb}$ (optional OCXO)
Graphical user interface	
Operating system	Windows 11 64-bit
Signals	
Signals	128 physical signals (inputs, outputs at 1 MS/s) 8 inputs and 8 outputs at 200MS/s 200 internal signals 24 simultaneous signals for data acquisition 8 simultaneous signals for 1 MS/s data streaming 4 simultaneous signals for 2 MS/s burst acquisition 24 simultaneous signals for real-time data logging
Data transfer	Via TCP/IP, 2kS/s default, up to 20kS/s, 1 MS/s × 8 channels for data streaming
Representation	32-bit floating point, real-world physical units

* Specifications are subject to change without notice.



- SPECSGROUP
- Sales and Service Partner

SPECS Surface Nano Analysis GmbH
Voltastrasse 5
13355 Berlin / Germany
www.specs-group.com

T +49 30 46 78 24-0
F +49 30 46 42 083
E info@specs-group.com

Version 02.26

SPECSGROUP