

NanoESCA • MARIS

HIGH PERFORMANCE REAL AND MOMENTUM SPACE ANALYSIS OF SMALL SAMPLES

KEY FEATURES

- Fully modular instrument
- Lateral resolution < 35 nm
- Momentum space resolution better than $5 \text{ m}\text{\AA}^{-1}$
- PEEM, MM and ARPES operation
- Imaging spin filter optionally



SPECS Surface Nano Analysis GmbH

T +49 30 46 78 24-0
E info@specs-group.com
H www.specs-group.com

 **Focus**™

A member of SPECSGROUP

